**REMSTIM<sup>®</sup>**  
RAPID EYE MOVEMENT STIMULATOR  
4000

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# Operating Instructions for the REMSTIM 4000

(Status 21.04.2020: subject to changes)

## Contents

Customer Service.....	23
Before the first use .....	23
Important Information .....	23
Scope of delivery.....	24
Brief instructions: Building up & Positioning of the REMSTIM 4000 .....	25
Setting up the REMSTIM 4000 .....	28
Starting the REMSTIM 4000 .....	29
Operation of the stimulation unit REMSTIM 4000.....	30
Functions.....	30
Start and stop .....	31
Focus-1 / Focus-2.....	31
Switch duplicated pacing signal and clock signal backwards.....	32
Adjust speed.....	32
Reset speed to default value .....	33
Adjust brightness .....	33
Reset brightness to default value .....	34
Positioning of the REMSTIM 4000 .....	34
Power supply of the REMSTIM 4000 .....	37
Practice: the first session with yourself! .....	37
Coach yourself in six steps with the REMSTIM 4000.....	38
Maintain the REMSTIM 4000 in top condition.....	39
Warranty.....	39
There is warranty in accordance with the statutory provisions.....	39
Notice pursuant to German Battery Act (BattG).....	39
Replacing the batteries of the remote control .....	40
Notice regarding disposal of old appliances and electronic scrap .....	40
Electromagnetic fields .....	40

## Customer Service

Thomas Buhl Coaching Instrumente  
Am Dachsbau 59  
13503 Berlin  
Germany

Service-Hotline: 0049-(0)30-39 82 14 19  
E-Mail: info@remstim.com  
Internet: 4000.remstim.com  
Device-type: REMSTIM 4000

## Before the first use

The REMSTIM 4000 - Rapid Eye Movement stimulator was developed especially for the stimulation of rapid eye movements. The speed and brightness of the stimulation signal can be adjusted individually to the needs during a stimulation sequence. The eye displacement - how far the eyes look to the left and to the right during a stimulation sequence - is set by the selected distance between stimulation unit of the REMSTIM 4000 and the eyes of the client.

The REMSTIM 4000 has been specially designed to replace the usual manual waving in an EMDR session. Everything else remains the same: from preparation to the end of a session. Prepare your client for stimulation and make sure that your client only moves his eyes and keeps his head still during a stimulation sequence.

Please read the operating instructions carefully and completely before first use of the REMSTIM 4000 and store them well. When passing on the REMSTIM 4000, these instructions are to be handed over.

## Important Information

1. The REMSTIM 4000 is a device used for the stimulation of eye movements
2. The use of REMSTIM 4000 does not replace the work of a psychotherapist, physician or alternative medicine practitioner. Therefore, an ongoing treatment should not be interrupted or cancelled, and any treatment needed in the future should not be postponed or omitted entirely.
3. The application of the REMSTIM 4000 does not nullify medical directives. The responsibility lies entirely with the user.
4. The use of the REMSTIM 4000 with persons who suffer or have suffered particularly from the following diseases should be checked by the user for compatibility:
  - a. Epilepsy
  - b. Eye disorders,
  - c. Cerebro-organic diseases,
  - d. Addictive disorders,
  - e. Dissociative disorders,
  - f. Ego disorders.

5. It is the responsibility and decision of the user to begin, continue or terminate the application of the REMSTIM 4000.
6. The use of REMSTIM 4000 can temporarily cause intense emotions.
7. Use the REMSTIM 4000 only for the purpose for which it is made.
8. If the eye movements are found to be unpleasant, particularly during the first application, it requires the optimization of the personal device settings of displacement (correct distance) and / or speed and / or brightness. In rare cases of dizziness, nausea or eye pain, stop REMSTIM 4000 and adjust settings.
9. Do not expose REMSTIM 4000 to extreme weather conditions (rain, snow, sun, wind).
10. Protect the REMSTIM 3000 against water.

## Scope of delivery

Before setting up the REMSTIM 4000, please take a look at the individual parts:

1. Stimulation unit REMSTIM 4000
2. Base (stand). Be careful, the stand is heavy, it weighs 3 KG. This ensures that the REMSTIM 4000 is safe.
3. Stand rod, extendable
4. Gooseneck, flexible.
5. Clamp. Already mounted on the gooseneck
6. 2 cable clamps (narrow, 15 mm), already connected to the gooseneck
7. 2 cable clamps (wide, 20 mm), already attached to the tripod
8. Power plug (5V) with USB A output
9. Power cable (USB A to DC plug), 1m
10. Extension cable 3m
11. Powerbank (battery) with USB A - output incl. Charging cable
12. Infrared remote control incl. 2 batteries (AAA)
13. Manual

Setting up and starting the REMSTIM 4000 take only a few minutes. You do not need a tool!

## Brief instructions: Building up & Positioning of the REMSTIM 4000

	<p>Step 1:</p> <p>Prepare components.</p> <p>Put stand on the ground.</p>
	<p>Step 2:</p> <p>Screw the stand rod firmly to the large thread in the base.</p>
	<p>Step 3:</p> <p>Loosen the locking collar on the stand rod and extend the stand by approx. 30 cm by pulling out the inner rod. Fix the inner rod with the locking collar.</p>

	<p>Step 4:</p> <p>Gooseneck mount on the tripod. To do this, hold the pulled-out stand rod and screw it tight to the gooseneck. Then retract the inner stand rod, as the height is usually sufficient.</p>
	<p>Step 5:</p> <p>Turn the gooseneck 90 degrees. In this position, insert the stimulation unit of the REMSTIM 4000 into the clamp so that the tube protrudes only a few mm on the side of the socket for the power connection.</p>
 <p>Modus focus 2 with battery operation</p>	<p>Step 6a:</p> <p>Battery operation</p>

 <p>Preparation of mains operation</p>	<p>Step 6b: mains operation</p>
<p>Examples of positioning</p>	
 <p>REMSTIM 4000 in <b>stimulation position</b> (high)</p>	<p>REMSTIM 4000 in high stimulation position</p>
 <p>REMSTIM 4000 in <b>stimulation position</b> (middle)</p>	<p>REMSTIM 4000 in the middle stimulation position</p>

 <p>REMSTIM 4000 in stimulation position (low)</p>	<p>REMSTIM 4000 in low stimulation position</p>
 <p>REMSTIM 4000 in parking position</p>	<p>REMSTIM 4000 in parking position</p>

Please read the further details on setting up and positioning in the following text. There you will find important information and many tips for the correct use of the REMSTIM 4000.

## Setting up the REMSTIM 4000

Place the stand on the floor and screw the stand rod firmly to the large thread in the base. Hand-tight tightening is sufficient.

Mount the gooseneck on the tripod. To do this, loosen the hand-width locking ring on the stand rod and extend the stand by approx. 30 cm by pulling out the inner rod. Fix the tripod with the locking ring. Now screw the gooseneck to the tripod. To do this, hold the pulled-out stand rod firmly and then screw it tight to the gooseneck. Hand-tight, if firm tightening is sufficient.

The Connection of gooseneck and tripod should be kept in mind. During the session, the R4 will swing from its pacing position to its parked position and back. To do this, you usually turn the R4 back and forth at an angle of 90 degrees around the tripod axis.

If the gooseneck and the tripod are screwed tight, the o.g. Move as desired the inner stand rod. However, if the gooseneck and tripod are not tightly bolted, the gooseneck would turn on the tripod and come off over time. That must be avoided.

Bend the gooseneck 90 degrees so that the REMSTIM 4000 lettering on the clamp is horizontal. In this position, the stimulation unit of the REMSTIM 4000 can be inserted into the clamp.

The clamp has a wide and a narrow end to insert the stimulation unit. This allows the stimulation unit to be carefully inserted at the far end. Use the end of the stimulation unit with the power socket and the sticker. On the sticker there are 2 elongated markers for positioning. For proper alignment, one of the bar marks must be visible in the upper gap between the two jaws of the clamp, depending on whether the R4 is placed on the left or right hand of the client.

The end of the stimulation unit with the power socket must protrude only a few millimeters from the narrow end of the clamp. If correctly adjusted, the circulating hairline of the pacemaker sticker terminates with the clamp.

The clamp accurately holds the pacing unit in place. The shape of the clamp allows for a slight movement of the stimulation unit in the terminal when touched. This avoids damage to the stimulation unit.

Push the stimulation unit into the far end of the clamp with the end of the power socket and the sticker until the stimulation unit protrudes a few millimeters from the narrow end of the clamp.

## Starting the REMSTIM 4000

Plug the power adapter into the power outlet and connect the power cord to the stimulation unit. You can use the cable clamps on gooseneck and tripod to fix the cable. The cable clamps can be disconnected and clamped by hand.

As soon as the power supply has been established, the stimulation unit emits a red light on the back (on the side facing away from the client), signaling that it is ready for operation.

With the help of the infrared remote control, you can now operate the functions of the stimulation unit.

Surely you know an infrared remote control already from a TV. Since they are light signals that you send out with the IR remote control, it is important for the reliable implementation of the function commands to face the infrared remote control with the transmitting diode to the stimulation unit. The light in the infrared, which is invisible to our eye, must reach the stimulation unit.

The infrared remote control comes with 2 AAA batteries already in use. Change them regularly, as only sufficiently charged batteries provide reliable operation. It is also advisable to have spare batteries ready so that the batteries can be changed quickly even during a session if necessary.

## Operation of the stimulation unit REMSTIM 4000

The stimulation unit of the REMSTIM 4000 is only operated with the remote control. The remote control has the following key functions:

Key	Function
	Start and stop the stimulation sequence
	Focus-1: To refocus the attention of the client. Focus 2: For positioning
	Switching between the display of the full duplicate stimulation signal and the single middle clock-indicating LED.
	Increase and decrease the speed of the stimulation pulse
	Increase and decrease the brightness of the stimulation pulse
	Reset speed to the default value.
	Reset brightness to the default value.

## Functions

The keys you will most often use are the:

- Red button: Start and stop of the stimulation signal
- Star key: Focus 1 / Focus 2

That's because you usually start and stop many stimulation sequences during a session. Of course, you can also run the stimulation sequence for a full 90 seconds. However, this will rarely happen, as clients signal by (subliminal) signs before the end of the 90 seconds that a stimulation sequence can be ended.

Frequently, and especially with new clients who have little experience with rapid eye movements, you will often use the star key (focus 1) to refocus during a stimulation sequence. If the client is focused again, just restart the stimulation sequence.

A special feature of the REMSTIM 4000 is the initial acceleration of the stimulation signal at each start. Read more under "Adjusting speed".

## Start and stop

Press the red button to start the stimulation sequence. The middle LED signals the start of the stimulation sequence with a light of 1 second.

A stimulation sequence ends automatically after 90 seconds.

You can stop the stimulation sequence at any time by pressing the red button again. The stimulation signal then runs to the middle and ends with the illumination of the middle LED for 1 second.

Then you can immediately restart a stimulation sequence with the red button.

## Focus-1 / Focus-2

The star key is the second most important key and controls the focus functions 1 and 2. Both focus functions end automatically after 30 seconds, if you have not previously stopped them by pressing the red key, ie by starting a stimulation sequence, or actively be terminated.

- Press the star key once: Starts the function Focus-1.

Focus 1: The middle LED starts flashing fast. The attention of the client is thereby focused.

As soon as the client is unable to follow the stimulation signal with the eyes correctly during the stimulation sequence, press the star key once. Immediately the stimulation sequence is stopped and replaced by a fast flashing of the middle LED. The client thereby refocuses the signal. When the red button (Start) is pressed, the flashing lights turns into a steady light for the duration of 1 second. Hereafter a new stimulation sequence starts.

If the client is unable to follow the stimulation signal repeatedly, it is advisable to reduce the speed of the stimulation signal.

- Press the star key again during the operation of the Focus 1 function (2nd time): Starts the function Focus-2.

Focus 2: The center and the outermost LEDs are lit steadily. Thanks to these "position lights", you can now position the stimulation unit REMSTIM 4000 in front of the client's eyes. Please read more about positioning the REMSTIM 4000 under Positioning.

- Press the star key again while operating the Focus 2 function (3rd time): Exits the focus functions.

To exit the Focus 1 function, press the star key twice in succession. The Focus 2 function will still start, but will be stopped immediately by the second key press.

## Switch duplicated pacing signal and clock signal backwards

The AV key allows during a stimulation sequence to switch back and forth between the stimulation signal duplicated on the back and the reduced clock signal.

The client-facing front presents 59 green LEDs behind satined acrylic glass that emit the stimulation signal. On the back there are 13 green LEDs, which duplicate / reflect the stimulation signal. This makes it possible to control, during a stimulation sequence, whether the client can correctly follow the stimulation signal.

The backward stimulation signal may be reduced as needed to a clock indicating flashing (clock signal) of the middle LED.

Each stimulation sequence always starts with the back-duplicated stimulation signal.

The possibility of choosing between the duplicated stimulation signal and the reduced clock signal on the back was created because a completely mirrored stimulation signal, especially in experienced clients, is not always necessary and the rear light signal "can also be annoying"!

Nevertheless, there should always be a minimum of backward signal to control.

## Adjust speed

A special feature of the REMSTIM 4000 is the initial acceleration of the stimulation signal at each start. As with handwaving, the final speed of the stimulation signal is not abruptly executed from the beginning. Rather, the stimulation signal always starts at a slow speed (60 hz / m) and increases it within a few seconds to the target speed, which was automatically saved by the last application. This ensures that the client is always slowly and safely tuned and "picked up" from the beginning even with a very fast stimulation signal.

Use the +/- buttons to increase and decrease the speed of the stimulation signal.

- Each time the + button is pressed during stimulation, the current speed is immediately increased by 10 half cycles per minute until the highest speed of 180 half cycles per minute is reached.
- Each time the - button is pressed during stimulation, the current speed is immediately decreased by 10 half cycles per minute until the slowest speed of 60 half cycles per minute is reached.

Speed levels: The last used setting is saved.

Level	Half Cycles Per Minute	Comment
S1	5	The slow steps S1 to S6 are mainly used for special techniques that work with slow stimulus such as weaving, flash technique, absorption technique.
S2	10	
S3	20	
S4	30	
S5	40	
S6	50	
1	60	
2	70	
3	80	
4	90	
5	100	
6	110	
7	120	Default Value
8	130	
9	140	
10	150	
11	160	
12	170	
13	180	

#### Definition of speed:

- The speed is defined by the number of illumination of the two outer LEDs per minute.
- When the light signal passes from one outer LED 1 to the other outer LED 2 and back to LED 1, this is 1 cycle.
- When the light signal passes from one outer LED 1 to the other outer LED 2, this is 1 half cycle
- Thus, the speed is defined by the number of half-cycles per minute (hz / m)

The default speed is 120 half cycles per minute.

### Reset speed to default value

To reset the speed to the default value, press the blue button and the + button in sequence.

### Adjust brightness

Use the </> buttons to increase and decrease the brightness of the stimulation signal.

- Each time the < button is pressed during stimulation, the current brightness is immediately increased by 1 level until the highest brightness level 7 is reached.
- Each time the > button is pressed during stimulation, the current brightness is immediately decreased by 1 level until the lowest brightness level of 1 is reached.

Brightness levels: The last used setting is saved

Brightness levels	Comment
1	
2	
3	
4	Default value
5	
6	
7	

The brightness has been chosen so that at the brightness level 4 (standard) the stimulation signal is clearly visible in a moderately illuminated room.

Depending on the lighting conditions, please adjust the brightness level.

### Reset brightness to default value

To reset the brightness to the default value, press the blue button and the < button one after the other.

## Positioning of the REMSTIM 4000

The REMSTIM 4000 can be flexibly positioned in front of the client.

The following factors are decisive:

1. Position of the holding device: laterally near the left or right knee of the client
2. Height of the holding device: REMSTIM 4000 at eye level, or slightly higher or lower
3. Centered alignment of the REMSTIM 4000
4. Distance of the REMSTIM 4000 to the eyes of the client
5. Horizontal alignment of the REMSTIM 4000

Ad 1: Usually the coach and client sit on opposite sides each other during an EMDR session, laterally offset. Depending on which hand the coach waves with, the right or left knee of coach and client come close.

If the coach waves with his right hand, the REMSTIM 4000's tripod is placed on the left side of the client's left knee. Thanks to the gooseneck, the stimulation unit of the REMSTIM 4000 is positioned in front of the client.

If the coach waves with his left hand, the REMSTIM 4000's tripod will be placed on the right outside of the client's right knee. Thanks to the gooseneck, the stimulation unit of the REMSTIM 4000 is positioned in front of the client.

Please note here that "left" and "right" describes different sides, depending on whose position you look from. Because coach and client sit opposite each other.


Important: The stimulation unit of the REMSTIM 4000 must be rotated 180 degrees in the clamp around its own axis (Attention: Turn the stimulation unit in the clamp only in the direction of the two arrows, which you also find the sticker near the hairline, so the sticker will not be damaged). The side of the tube with the socket for the power connection remains in the terminal.

Why is a 180 degree turn required depending on left or right positioning?

If you look closely at the stimulation unit, you can see a large number of LEDs (59 LEDs) through the satin-finished acrylic glass of the tube on one side. This is the page that must always be facing the client to present a continuous stimulation signal.

On the back you will find 13 LEDs, which produce the duplicated stimulation signal for the coach, and the red LED, which indicates the operational readiness.

As a result, depending on the right or left hand positioning by the client, the stimulation tube must be rotated 180 degrees about its axis to align the 59 LEDs to the client.

However, you will typically need to do this only once, as I'm sure you will not change your session layout (left or right) from session to session.

Ad 2: The tripod allows you to adjust the height. However, you will usually leave the tripod at its lowest height as the gooseneck increases the stimulation unit.

The flexible gooseneck makes it possible to find an optimal position of the stimulation unit for the client. You can also bend the goose neck down like a side "S" to reduce the height.

It is up to you how high you position the stimulation unit in front of the client. The height has a significant influence on whether the client performs the rapid eye movements while he

- looks straight ahead
- looks up a bit
- looks down a bit

In addition, the height has an effect on the extent to which the direct eye contact with the client could be somewhat obscured by the stimulation unit. If the direct eye contact should be disturbed a little, a slight change in the height of the position of the stimulation unit immediately will help.

Tip: I recommend positioning the stimulation unit just above the client's eyes, as this ensures that the

eye is completely undisturbed while also providing a clear view of the client's face; every change of facial expression of the client is thus clearly visible.

Ad 3: The stimulation unit should be in a centered position in front of the client's eyes. Activate the Focus-2 function for centering by pressing the star key on the remote control twice. The REMSTIM 4000 must be connected to a power source.

The focus 2 function ends after 30 seconds and can then be restarted immediately by pressing the star key twice. (see also under Functions / Star key)

The center and the outermost LEDs are now lit on the stimulation unit. While the client is sitting in his usual sitting position, position the stimulation unit with the center LED in the center of the client's field of view.

to 4: For an optimal distance between the stimulation unit and the eyes of the client, the stimulation unit should be located about 30 centimeters (that is the length of a commercial ruler) in front of the eyes (maybe slightly higher, so that the client looks a little bit upwards). Now ask the client if the distance is comfortable for him and if he can see the glowing LED sharply.

If the LED appears too close to the client, or if he / she is unable to see it clearly, increase the distance between the stimulation unit and the client's eyes in small increments while constantly asking for visual acuity. To do this, move the entire holding device together with the stand and stimulation unit.

The distance between the stimulation unit and the client's eyes should not be too large, because even if the REMSTIM 4000 offers a generous stimulation range of 58 cm between the outermost LEDs, a long distance between the stimulation unit and the client's eyes may cause an insufficient movement of the eyes.

Test the eye moving range: Then ask the client to focus alternately on the left and then the right side of the illuminated LED without moving the head. This will test whether the range of the eyes (that is, how far the eyes move left and right during a stimulation sequence) is OK. If the range is uncomfortable, increase the distance between the stimulation unit and the client's eyes in small increments while constantly asking for appropriateness.

Ad 5: Make sure that the stimulation unit is aligned horizontally in front of the client. In the optimal case all LEDs are at the same height, the stimulation unit "does not hang diagonally" in front of the client. You can check this well thanks to the Focus-2 function.

Note: the stimulation unit clamp has a pivot that is fixed with a locking screw. It can be loosened with a screwdriver and tightened.

I recommend having this screw tightened. All necessary adjustments regarding the horizontal orientation of the stimulation unit can be made very precisely with the flexible gooseneck. To do so, hold the short shaft of the clamp just below the hinge and a piece of the upper part of the gooseneck with your hand. With the other hand, hold the lower part of the gooseneck, which is screwed to the tripod. Now bend the gooseneck upon need and align the stimulation unit horizontally.

## Power supply of the REMSTIM 4000

The REMSTIM 4000 is operated with 5V DC. Please use only the supplied components for power supply. There are 2 options available.

1. Mains operation: Mains plug with USB A connection
2. Battery operation: Powerbank (battery incl. Charging cable), allows about 40 hours of continuous operation with fully charged battery.

Disconnect the REMSTIM 4000 from the power supply while it is not in use.

To 1: To use the power plug, connect the supplied adapter cable (USB A to DC plug, 1m) to the power adapter. Then extend the adapter cable with the extension cable (3m). Then connect the extension cable to the power socket of the REMSTIM 4000. As soon as you insert the power plug into a socket, the REMSTIM 4000 is ready for use, which is signaled by the red LED on the back of the stimulation unit lighting up.

The cable can be guided along the tripod and gooseneck with the cable clamps.

Ad 2: To operate the REMSTIM 4000 with the power bank (battery), connect the power bank with the supplied adapter cable (USB A to DC plug, 1m). Then connect the DC plug of the adapter cable to the power socket of the REMSTIM 4000. The REMSTIM 4000 is now ready for use, which is signaled by the red LED on the back of the stimulation unit lighting up.

The cable can be guided along the tripod and gooseneck with the cable clamps.

The power bank can be attached to a cable clamp (eg at the bottom) thanks to the key chain ring. Loosen the clamp to fix the ring as in the picture.

The supplied power bank (battery) has been partially charged for transport. For test purposes, the charge of the battery is sufficient. Before the first regular operation of the REMSTIM 4000 it is necessary to fully charge the battery with the supplied charging cable.

## Practice: the first session with yourself!

I recommend that you first use the REMSTIM 4000 in an "EMDR self-coaching session" with yourself once you set it up and put it into operation.

Sit on the client chair. In this way, you will directly learn how to best position the REMSTIM 4000 in front of your eyes. The following self-coaching session gives you the opportunity to practice the handling and operation of the REMSTIM 4000 with the infrared remote control.

You will notice that the red key (Start / Stop) and the star key for the focus functions (especially focus 1) are the most commonly used keys.

## Coach yourself in six steps with the REMSTIM 4000

Please note that the EMDR self-coaching process differs from the EMDR standard protocol. Learn more about self-coaching with EMDR at [remstim.com](http://remstim.com).

Here is a quick start guide that will allow you to use the REMSTIM 4000 with yourself for the first time.

### 1. Step: Your subject

- What moves you the most at this moment?
- What kind of thoughts are burdening you currently? What is the most unpleasant image, the most unpleasant sentence?
- What's on your mind?
- Concentrate on the idea, the situation, the thing, the image, the phrase, etc. that burdens you the most!

### 2. Step: Your emotion

- Do you feel something while you concentrate on your subject?
- Which emotion do you feel?
- Is it for example: Anger, resentment, outrage, annoyance, fear, anxiety, grief, loss, helplessness, powerlessness, surprise, confusion, "being in the wrong movie", guilt, responsibility, shame, humiliation, insult, degradation, disgust, hatred, abhorrence, repugnance, "horror", fatigue, thirst, hunger, inability to feel? Or another emotion? (The list represents a number of unpleasant emotions and does not claim to be complete)

### 3. Step: Body scan I – Your physical sensations

- Do you feel something in your body? Is there a body echo of your emotion?
- Do you feel an emotion, a physical sensation somewhere in a part of your body?
- Can you perceive some kind of pressure, drag, heat, cooling or the like within you?
- Or do you perceive nothing, and it is like a non-sensation?

### 4. Step: The intervention (Processing)

- Feel your topic and your emotion(s) within you. Feel your physical sensation (including non-sensation) now!
- Now start the application of REMSTIM 3000
- Breathe gently and regularly!
- When the REMSTIM 3000 application ends, take a deep breath and let your feelings pour out of you when you exhale.
- What has changed internally? Observe your emotion(s), your body feeling, when you now think about your topic from the beginning (associations?).

→ Go back to the first point of intervention (4th step) and begin the process again. Go through the process as often as you want. For example, until you have distanced yourself from your topic and feel liberated.

### 5. Step: Body scan II

- Think about your subject. Sense and search for small residual feelings. If you find any, intervene with the REMSTIM 3000 as described in the 4th step until the best possible result is achieved.

### 6. Step: Finally – The anchoring

- Now think of your topic and notice the first positive thought that comes to your mind.
- Repeat this for your pleasant thoughts in the mind again and again. Intervene now with the REMSTIM 3000 – Choose the slow chaser light speed.
- Breathe gently and regularly and repeat your pleasant thought.

After using the REMSTIM 3000, breathe deeply and release your feelings as you exhale.

## Maintain the REMSTIM 4000 in top condition

- Clean the surface of REMSTIM 3000 regularly. To do this, only use a dry cloth (glasses cleaning cloth).
- Damage to the device: If the device is broken or damaged, do not attempt to repair the device yourself. Contact the customer service department in case of damage (see Customer Service section).

## Warranty

There is warranty in accordance with the statutory provisions.

## Notice pursuant to German Battery Act (BattG)

The infrared remote control of the REMSTIM 4000 is powered by 2 batteries, which are included in the delivery. Therefore, we are obliged under the German Battery Act (BattG) to draw your attention to the following:

Batteries and accumulators must not be disposed of with household waste. You are rather legally obliged to return used batteries and accumulators!

Old batteries can contain harmful substances that could harm the environment or your health if improperly stored or disposed of. Batteries also contain important raw materials such as iron, zinc, manganese, or nickel, and can be recycled.

After using the battery, you can either

- return it to us (Thomas Buhl Coaching Instrumente, Am Dachsbau 59, 13503 Berlin, Germany)
- or return it in your area free of charge (for example, in a retail market or at communal collecting points).

The return of batteries to retail outlets is limited to end users for the disposal of usual quantities, and of those old batteries which the retailer offers or has previously offered in its product range.

The cross-out waste bin sign means that you cannot dispose of batteries and accumulators with household waste. The signs under the waste bin on batteries stand for:

- Pb: Battery contains lead
- Cd: Battery contains cadmium
- Hg: Battery contains mercury

## Replacing the batteries of the remote control

Insufficient battery power will cause poor infrared transmit power, causing unreliable remote control commands to the REMSTIM 4000 stimulation unit. If the performance of the supplied batteries (2 x LR03 AAA) wears off, replace the old batteries with new ones.

To do this, loosen the cover of the battery compartment on the back of the remote control by sliding it downwards. Please do not use a wrench or similar tool.

Remove the old batteries and insert the new batteries. Close the battery compartment.

## Notice regarding disposal of old appliances and electronic scrap

From March 2006 onwards, old devices may no longer be disposed of with household garbage as per the Electrical and Electronic Equipment Act (ElektroG). This not only includes the typical electrical appliances such as televisions, washing machines or refrigerators, but also video games, speakers, computers, laptops, MP3 players, peripherals, printers, telephones, electric-powered children's toys or other electric-powered old devices. These devices can be recognised through the symbol of a "crossed-out garbage can". The symbol can be found on the equipment or packaging.

The electronic scrap can be given free of cost to municipal collection points. Your municipal waste management or the consumer centres can give you information on the exact location of the collection points.

## Electromagnetic fields

The EMDR device REMSTIM 4000 complies with all relevant standards and regulations regarding exposure to electromagnetic fields.

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<p>REMSTIM 4000 in <b>Stimmulationsposition</b> (tief)</p> <p>The image shows the REMSTIM 4000 device in its deep stimulation position. The chair is tilted back, and the footrest is extended. A metal stand with a horizontal bar is positioned behind the chair, and a small white table is to the right.</p>	<p>REMSTIM 4000 en position de stimulation profonde.</p>
<p>REMSTIM 4000 in <b>Parkposition</b></p> <p>The image shows the REMSTIM 4000 device in its parking position. The chair is upright, and the footrest is retracted. The metal stand and table are in the same position as in the previous image.</p>	<p>REMSTIM 4000 en position de parking.</p>

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